

SAFETY MED LTD

Unit 1A, Hallam Business Centre, Hallam Street, Heaviley, Stockport,
Cheshire. SK2 6PT. Tel: 01614773163. Fax: 01614774380.

E-mail: sales@safetymedltd.co.uk

www.safetymedical.co.uk



INDEX

ISOLATION UNIT	2
MOBILE OR STATIC EXTENSION UNITS	3
IV EXTENSION UNITS	4, 5
SURGEONS STOOLS	6
INSTRUMENT AND DRESSING TROLLEYS	7
TRANSFUSION STANDS	8
BOWL STANDS	9
SOILED LINEN TROLLEY	9
MAYO TABLES	10
PLATFORM STEPS	11
SURGEONS STEPS	12
BARIATRIC HEAVY DUTY STEPS	13
CLAMPS	14
ARM TABLE AND X RAY TABLES	15
ANAESTHETIC POLES	16
PELVIC SUPPORTS	16
BUTTOCK SUPPORTS	16
LATERAL SUPPORTS	17
CARTER BRAIN SUPPORT	17
LITHOTOMY POLES CRANKED AND CURVED	18
LLOYD DAVIES SUPPORTS	19
TABLE SIDE RAILS	19
PATIENT RESTRAINT STRAPS WITH CLAMPS	19
ADULT KNEE CRUTCH SUPPORTS	20
CHILD KNEE CRUTCH SUPPORTS	20
BARIATRIC SUPPORT UNIVERSAL	21
LATERAL SUPPORT (SPECIAL)	22
LEG-PREP SUPPORT	23
HEPWORTH SUPPORTS	24
ARM RETAINERS	25
STAINLESS STEEL SHELVING UNITS	26
PLASTER TROLLEY	27
SPECIALS (TAILOR MADE)	29-34
MATTRESS AND UPHOLSTERY	35- 45

- WE MARKET A WIDE RANGE OF HOSPITAL THEATRE EQUIPMENT AND ACCESSORIES.
- WE ALSO SUPPLY A VARIED RANGE OF HOSPITAL EQUIPMENT FOR USE IN THE GENERAL WARDS.
- DUE TO OUR IN HOUSE DESIGN FACILITIES WE CAN PROVIDE A COMPREHENSIVE SERVICE TO SUPPLY TAILOR MADE EQUIPMENT TO MEET YOUR SPECIFICATION.
- WE HAVE EXTENSIVE EXPERIENCE OF WORKING WITH MEDICAL PROFESSIONALS AND CAN THEREFORE PROVIDE ADVICE AND SUGGESTIONS IN ORDER TO ACHIEVE YOUR REQUIREMENTS.
- ALL ITEMS ARE PRODUCED USING MATERIALS TO THE HIGHEST STANDARD.
- ALL UPHOLSTERY IS MANUFACTURED FROM TOP QUALITY ANTISTATIC AND ANTIBACTERIAL FABRIC.

UNIQUE MOBILE SAFETY ISOLATION EXTENSION UNIT



THE ESSENTIAL SAFETY DEVICE

BRINGING THE POWER TO THE TABLE, MOBILE EXTENSION UNIT WITH 8 X 13 AMP SWITCHED SOCKETS, POWERED BY 2.9KVA SAFETY ISOLATION TRANSFORMER. PATIENT & USER FRIENDLY. POWER INDICATOR & OVERLOAD RE-SET CIRCUIT BREAKER VISIBLE VIA CLEAR PETG TOP COVER, WHICH PROTECTS SOCKETS FROM WET CONDITIONS. SUPPLIED WITH 5M LEAD AS STANDARD. ALL STAINLESS STEEL CONSTRUCTION (SATIN FINISH). ANTISTATIC CASTORS. FULLY CE MARKED & BUILT TO IEC6060-1 STANDARD.

OPERATING THEATRE & INTENSIVE CARE

MOBILE OR STATIC EXTENSION UNITS

EXTENSION UNIT
8 X 13 AMP SOCKETS &
5M LEAD WITH
A/S CASTORS OR FEET.



EXTENSION UNIT
10 X 13 AMP SOCKETS &
5M LEAD WITH
A/S CASTORS OR FEET.

LOW LEVEL EXTENSION UNIT
8 X 13 AMP SOCKETS &
5M LEAD WITH
A/S CASTORS OR FEET.
CLEAR PETG TOP COVER.



MANUFACTURED IN 304 STAINLESS STEEL

COPYRIGHT ©

OPERATING THEATRE & INTENSIVE CARE

BRING THE POWER TO THE TABLE



5 LEG (HEAVY) IV STAND WITH 5 SIDED, 10 WAY
X 13 AMP SOCKETS & 5M LEAD

4 LEG (HEAVY) IV STAND WITH 4 SIDED, 8 WAY
X 13 AMP SOCKETS & 5M LEAD

USE AS CHARGING STAND FOR
8 TO 10 PUMPS ETC.

MOBILE OR STATIC EXTENSION UNIT
WITH 10 X 13 AMP SOCKETS &
OR
8 X 13 AMP SOCKETS
& 5M LEAD



MANUFACTURED IN 304 STAINLESS STEEL

COPYRIGHT ©

USE AS CHARGING STAND FOR PUMPS ETC.



**SPECIAL SUPPLIED TO SOLVE PROBLEM AT
UNIVERSITY HOSPITAL COVENTRY AND WARWICK.
DESIGNED WITH HELP FROM MEBS DEPT.
TRANSFUSION STAND WITH LOW CENTRE OF GRAVITY,
EXTRA WEIGHT & 10 WAY 13 amp EXTENSION UNIT LOCATED
OVER IV POLE, ALLOWING UPTO 6 PUMPS TO BE CLAMPED AT BED
MATTRESS HEIGHT, WITH COMMON EARTH TAG TO EXTENSION UNIT AND IV STAND.**

MANUFACTURED IN 304 STAINLESS STEEL

COPYRIGHT ©

OPERATING THEATRE ACCESSORIES

SURGEONS TUBULAR STOOLS & CHAIRS

ALL UPHOLSTERY MANUFACTURED FROM TOP QUALITY ANTISTATIC & ANTIBACTERIAL FABRIC



FIVE LEG STOOL WITH
A/S CASTORS
OR FEET



FOUR LEG STOOL WITH
A/S CASTORS
OR FEET



FIVE LEG STOOL & BACKREST
WITH
A/S CASTORS
OR FEET



FOUR LEG STOOL & BACKREST
WITH
A/S CASTORS
OR FEET

REPLACEMENT STOOL COVER 13" dia. REPLACEMENT BACKREST PAD
BRAKED CASTORS.

MANUFACTURED IN 304 STAINLESS STEEL.

OPERATING THEATRE ACCESSORIES

INSTRUMENT & DRESSING TROLLEYS



FLUSH WELDED SHELVES
& A/S CASTORS
18" X 18"
18" X 24"
18" X 30"
18" X 36"



3 LIP WELDED SHELVES
& A/S CASTORS
18" X 18"
18" X 24"
18" X 30"
18" X 36"



FLAT LOOSE SHELVES, REVERSIBLE
OPTIONAL GUARD RAILS.
A/S CASTORS
18" X 18"
18" X 24"
18" X 30"
18" X 36"



3 LIP LOOSE SHELVES, REVERSIBLE
OPTIONAL GUARD RAILS.
A/S CASTORS
18" X 18"
18" X 24"
18" X 30"
18" X 36"

OPERATING THEATRE ACCESSORIES

TRANSFUSION STANDS

FOUR LEG TRANSFUSION STAND WITH
2 HOOKS & A/S CASTORS.

LIGHT OR
HEAVY

EXTRAS.

4 HOOKS
4 DOUBLE HOOKS
4 PIGTAIL HOOKS



FIVE LEG TRANSFUSION STAND WITH
2 HOOKS & A/S CASTORS.

LIGHT OR
HEAVY

EXTRAS.

4 HOOKS
4 DOUBLE HOOKS
4 PIGTAIL HOOKS

HEIGHT 53" MIN. TO 84" MAX.

OPERATING THEATRE ACCESSORIES

BOWL STAND (single)
BOWL STAND (double)
NOT ILLUSTRATED
SOILED DRESSING CONTAINER
STAINLESS STEEL BOWL 14" dia
BUCKET & STAND
STAINLESS STEEL BUCKET



OVERBED BOWL STAND (single)
OVERBED BOWL STAND (double)

SOILED LINEN TROLLEY (single)
SOILED LINEN TROLLEY (double)
NOT ILLUSTRATED
SOILED LINEN TROLLEY (with platform)



OPERATING THEATRE ACCESSORIES



MAYO TABLE
TRAY SIZE 310mm x 520mm
12" x 21-5"
STAINLESS STEEL



MAYO TABLE HYDRAULIC
TRAY 540mm x 740mm
21" x 29"



MAYO TABLE HYDRAULIC
TABLE 30° LEFT & RIGHT



MAYO TABLE HYDRAULIC
TABLE HEIGHT TRAVEL 190 – 200mm
FOOT PUMP ACTION

OPERATING THEATRE ACCESSORIES



PLATFORM STEPS
ANTI SLIP FEET
ALL STAINLESS STEEL
SIZE 300mm x 460mm x 130mm HIGH
12" x 18" x 5"

SUPPORTING RIBS ON
UNDERSIDE



PUSH TOGETHER TO FORM A
LARGER PLATFORM

OPERATING THEATRE ACCESSORIES



SURGEONS STEPS

18 x 12 x 6"

30 x 12 x 6"

30 x 12 x 8"




OPERATING THEATRE ACCESSORIES




HEAVY DUTY STEPS
(BARIATRIC)




OPERATING THEATRE TABLE ACCESSORIES




ROTARY TOOTH CLAMP TO FIT
16mm OR 5/8" BAR
UK CLAMP
EUROPEAN CLAMP



UNIVERSAL CLAMP FOR ROUND
OR FLAT BAR
UK CLAMP
EUROPEAN CLAMP



END ON ROUND CLAMP
16mm OR 5/8" BAR
UK CLAMP
EUROPEAN CLAMP



END ON FLAT CLAMP
UK CLAMP
EUROPEAN CLAMP

OPERATING THEATRE TABLE ACCESSORIES

NARROW ARM TABLE 24" X 8"
WITH 16mm dia PIN

ANTISTATIC PAD



SQUARE ARM TABLE 17" X 17"
WITH 2 X 16mm dia PINS

ANTISTATIC PAD

LARGE X-RAY ARM TABLE 27" X 19"
WITH 2 X 16mm dia PINS & TELESCOPIC LEG

ANTISTATIC PAD



ALL UPHOLSTERY MANUFACTURED FROM TOP QUALITY ANTISTATIC & ANTIBACTERIAL FABRIC.
MANUFACTURED IN 304 STAINLESS STEEL. ALL PRICES EX WORKS PLUS VAT

OPERATING THEATRE TABLE ACCESSORIES

ANAESTHETIC SCREEN POLE
METALWORK ONLY
ANAESTHETIC SCREEN
PAD



PELVIC SUPPORT
METALWORK ONLY
PELVIC SUPPORT PAD

BUTTOCK SUPPORT
METALWORK ONLY
BUTTOCK SUPPORT PAD



ALL UPHOLSTERY MANUFACTURED FROM TOP QUALITY ANTISTATIC & ANTIBACTERIAL FABRIC.
MANUFACTURED IN 304 STAINLESS STEEL. ALL PRICES EX WORKS PLUS VAT.

OPERATING THEATRE TABLE ACCESSORIES



LATERAL SUPPORT CURVED
METALWORK ONLY
LATERAL SUPPORT PAD

LATERAL SUPPORT (TUDOR EDWARDS)
METALWORK ONLY
LATERAL SUPPORT PAD



ARM SUPPORT (CARTER BRAINE)
METALWORK
ARM SUPPORT PAD

ALL UPHOLSTERY MANUFACTURED FROM TOP QUALITY ANTISTATIC & ANTIBACTERIAL FABRIC.
MANUFACTURED IN 304 STAINLESS STEEL. ALL PRICES EX WORKS PLUS VAT.

OPERATING THEATRE TABLE ACCESSORIES



LITHOTOMY POLES CRANKED (pair)
METALWORK ONLY
SLEEVE 18" LONG
SLEEVE 22" LONG
STRAPS PADDED WITH
BILLIT CLIPS (pair)



LITHOTOMY POLES CURVED (pair)
METALWORK ONLY
SLEEVE 18" LONG
SLEEVE 22" LONG
STRAPS PADDED WITH
BILLIT CLIPS (pair)

OPERATING THEATRE TABLE ACCESSORIES

LLOYD DAVIES LEG SUPPORTS
METALWORK ONLY
LLOYD DAVIS PADS (pair)
LLOYD DAVIS LEG SUPPORT
STRAPS



OPERATING TABLE SIDE RAIL
METALWORK ONLY
SIDE RAIL PAD

PATIENT RESTRAINT STRAPS (pair)
WITH CLAMPS
CLAMPS ONLY (pair)



ALL UPHOLSTERY MANUFACTURED FROM TOP QUALITY ANTISTATIC & ANTIBACTERIAL FABRIC.
MANUFACTURED IN 304 STAINLESS STEEL. ALL PRICES EX WORKS PLUS VAT.

OPERATING THEATRE TABLE ACCESSORIES



ADULT KNEE CRUTCH SUPPORTS

PAIR COMPLETE
PADS (PAIR)



CHILD KNEE SUPPORTS
PADS (PAIR)



MANUFACTURED IN 304 STAINLESS STEEL. ALL PRICES EX WORKS PLUS VAT.
ALL UPHOLSTERY MANUFACTURED FROM TOP QUALITY ANTISTATIC & ANTIBACTERIAL MATERIAL

OPERATING THEATRE TABLE ACCESSORIES




BARIATRIC SUPPORT UNIVERSAL
CLAMPS NOT INCLUDED



16mm PINS STANDARD MODEL
ROLL 600mm LONG x 100mm DIA

SUPER BARIATRIC SUPPORT (NOT SHOWN)
PINS 16mm or 18mm
ROLL 600mm LONG x 150mm DIA
HEAVY DUTY CENTRE SHAFT

OPERATING THEATRE TABLE ACCESSORIES



LATERAL SUPPORT
(SPECIAL)
COMPLETE

CLAMPS AND CABLES
NOT INCLUDED



PADS SET
OF THREE

MANUFACTURED IN 304 STAINLESS STEEL. ALL PRICES EX WORKS PLUS VAT.
ALL UPHOLSTERY MANUFACTURED FROM TOP QUALITY ANTISTATIC & ANTIBACTERIAL MATERIAL.

OPERATING THEATRE ACCESSORIES



LEG PREP SUPPORT
5 degree increments



ADJUSTABLE ANKLE
WIDTH AND DEPTH



ANTISTATIC AND ANTI BACTERIAL PADS
FOR COMFORT

MANUFACTURED IN 304 STAINLESS STEEL. ALL PRICES EX WORKS PLUS VAT.
ALL UPHOLSTERY MANUFACTURED FROM TOP QUALITY ANTISTATIC & ANTIBACTERIAL MATERIAL.

OPERATING THEATRE ACCESSORIES



HEPWORTH SUPPORTS
X RAY TRANSLUCENT
POLYCARB/PET
MATERIAL

STANDARD UNIT
PAIR



SMALL UNIT
PAIR

LARGE BASE UNIT (NOT SHOWN) PAIR

POLYCARBONATE ARM & LEG RETAINERS & ANAESTHETIC TUBE HOLDERS



L TYPE 5" HIGH ARM RETAINER
L TYPE 8" HIGH ARM RETAINER
J TYPE ARM RETAINER
STEP TYPE ARM RETAINER
FLAT ANAESTHETIC TUBE HOLDER
L TYPE ANAESTHETIC TUBE HOLDER
(not illustrated)
ALL ARM RETAINERS 4" WIDE
ANAESTHETIC TUBE HOLDERS 6" WIDE



TEA-COSY PAD 5" ARM RETAINER
TEA-COSY PAD 8" ARM RETAINER

ALL UPHOLSTERY MANUFACTURED FROM TOP QUALITY ANTISTATIC & ANTIBACTERIAL MATERIAL.
ALL PRICES EX WORKS PLUS VAT.

STAINLESS STEEL SHELVING UNITS



6 SHELF UNIT
LARGE CASTORS – 2 BRAKED
SIZE 446mm x 910mm x 1520mm HIGH
SIZE 18" X 36" X 60" HIGH
EASY TO CLEAN

TAPERED SHELF UNIT WITH 4 SHELVES
AND ONE TOP SHELF
SIZES AS ABOVE



EASY TO CLEAN
(NO HIDDEN CORNERS; WHOLE UNIT)



STAINLESS STEEL PLASTER TROLLEY'S



PLASTER TROLLEY
TWO DRAWER
460mm x 900mm x 860mm HIGH
18" X 36" X 34" HIGH
TWO BRAKED CASTORS

STAINLESS STEEL BUCKET ON
SWIVEL BRACKET



SPECIALS



WORK STATION – DIVIDER FOR RECOVERY UNIT

SPECIALS



PUMP CHARGING STAND
12 x 13 AMP SOCKETS
12 PUMP POSITIONS (WALL FIXING)



MANUFACTURED IN 304 STAINLESS STEEL. ALL PRICES EX WORKS PLUS VAT

SPECIALS



MODIFICATION FOR OPERATING TABLE HEAD END



SPECIALS



MODIFICATION FOR OPERATING TABLE HEAD END



SPECIALS



ROTARY TOOTH CLAMPS (LARGE)
TO TAKE 1" DIAMETER CAMERA POLE

SPECIALS

EXTENSION UNIT 10 WAY
MOUNTED ON 5 LEG HEAVY DUTY STAND
(WAIST HEIGHT)



SPECIALS



MASK BOX HOLDER
(WALL MOUNTED)



GLOVE BOX HOLDER
(WALL MOUNTED)



OPERATING TABLE MATTRESSES

<u>CODE.</u>	<u>DESCRIPTION.</u>
SM001	Mattress Compatible with Model M 3-section mattress; comprising of hinged trunk section with perineal cut out, head & leg sections. 75 mm Deep.
SM001M	Mattress with 25mm memory foam topper
SM002	Mattress Compatible with Model MR 3- section mattress; comprising of hinged Trunk section with perineal cut out, head & leg sections. 75 mm Deep.
SM002M	Mattress with 25mm memory foam topper
SM003	Mattress Compatible with Model MS-M+ 3- section mattress; comprising of hinged trunk section with perineal cut out, head & leg sections. 75 mm Deep.
SM003M	Mattress with 25mm memory foam topper
SM004	Mattress Compatible with Models J & K. 3- section mattress; comprising of main trunk section with perineal cut out, head & leg sections. 75 mm Deep
SM004M	Mattress with 25mm memory foam topper
SM005	Mattress Compatible with Model J3 4- section mattress; comprising of upper trunk with perineal cut out, lower trunks with perineal fill in, head & leg sections 75 mm Deep.
SM005M	Mattress with 25mm memory foam topper

OPERATING TABLE MATTRESSES

<u>CODE.</u>	<u>DESCRIPTION.</u>
SM006	Mattress Compatible with Models GX & GXL 3- section mattress; comprising of head & upper Trunk with perineal cut out, hinged lower trunk & leg section with perineal fill in. 75 mm Deep.
SM006M	Mattress with 25mm memory foam topper
SM007	Mattress Compatible with Belfast Tables 3- section mattress; comprising of trunk section with perineal cut out, lower trunk & head section. 75 mm Deep.
SM007M	Mattress with 25mm memory foam topper
SM008	Mattress Compatible with Model X3 3- section mattress; comprising of trunk section with perineal cut out, head & leg sections. 75 mm Deep.
SM008M	Mattress with 25mm memory foam topper
SM009	Mattress Compatible with model MXPT. 5- section mattress; comprising of hinged upper & lower trunk, head & leg sections. 75 mm Deep.
SM009M	Mattress with 25mm memory foam topper
SM010	Mattress Compatible with International 5 5- section mattress; comprising of hinged upper & lower trunk, head & leg sections. 75 mm Deep.
SM010M	Mattress with 25mm memory foam topper
SM011	Mattress Compatible with Stille 3000 5 – section mattress; comprising of upper /middle lower trunk section with Perineal cut out, head & leg sections. 75 mm Deep.
SM011M	Mattress with 25mm memory foam topper

OBSTETRIC MATTRESSES

<u>CODE.</u>	<u>DESCRIPTION.</u>
SM101	Mattress Compatible with Rocket MK I 2- section mattress; comprising of upper section A & lower section B.
SM101M	Mattress with 25mm memory foam topper
SM102	Mattress Compatible with Rocket MK II 2- Section mattress; comprising of upper hinged section A & lower section B.
SM102M	Mattress with 25mm memory foam topper
SM103	Mattress Compatible with Birthday Bed 3- section mattress; comprising of upper trunk section A, middle trunk section B & lower trunk section C.
SM103M	Mattress with 25mm memory foam topper
SM104	Mattress Compatible with Birthright Bed. 2- section mattress comprising of upper trunk & lower trunk.
SM104M	Mattress with 25mm memory foam topper

PATIENT TROLLEY MATTRESSES

SM201	Patient Trolley mattress flat, with no hinge can be manufactured to any size required by the customer e.g:183 cm x 56 cm x	50 mm	75 mm
SM201M	Mattress with 25mm memory foam topper	50mm	75mm
SM202	Patient Trolley mattress with hinge can be manufactured to any size required by the customer. E.g.: 183 cm x 56 cm	50 mm	75 mm
SM202M	Mattress with 25mm memory foam topper	50mm	75mm

COT SIDE COVERS

Cot side covers can be manufactured to any size or design, using a selection of materials. Specification sheets and material samples are available on request, to aid in customer ordering.

SMALL

To include sizes 1070mm - 1270mm length

Single Sided with Velcro Flaps

pair

Tea Cosy Type

pair

MEDIUM

To include sizes 1270mm - 1530mm length

Single Sided with Velcro Flaps

pair

Tea Cosy Type *pair*

LARGE

To include sizes 1530mm onwards

Single Sided with Velcro Flaps

pair

Tea Cosy Type

pair

OPERATING THEATRE ACCESSORIES

CODE.

DESCRIPTION.



SM301

Rubens Type Pillow

Adult & child size.



SM302

Cardiff Head Stabiliser



SM303

Whitlock Head Rest



SM304

Horse Shoe Head Rest

SM305

Anaesthetic Head Support

SM306

Anti-Static Covered Pillow

OPERATING THEATRE ACCESSORIES



SM307 **Doughnut Head Ring - (Adult)**

SM308 **Doughnut Head Ring - (Child)**

SM309 **Tea Cosy Pad for Arm Retainer**
Please state for size required.



SM310 **Arm Board Pad**
Please state size required.



SM311 **Guttered Arm Board Pad**
Please state size required.



SM312 **Turret Heel Support Pair**



SM313 **Turret Heel Support Double**

OPERATING THEATRE ACCESSORIES



SM314

D Shaped Pillow

305mm x 100mm x 65mm deep

SM315

D Shaped Pillow

308mm x 100mm x 65mm deep

SM316

D Shaped Pillow

480mm x 100mm x 65mm deep

SM317

Seat Pad for Anaesthetic Chair

Please state size required

SM318

Back Rest Pad for Anaesthetic Chair

Please state size required

SM319

Stool Cushions

305mm dia; 330mm dia; 355mm dia.

OPERATING THEATRE ACCESSORIES



SM320

Lloyd Davies
Crutch Pads pair.



SM321

Buttock Support
Tubular Pad 255mm x 100mm dia.

SM322

Pelvic Support
Tubular Pad 280mm x 100mm dia.



SM323

Lateral Support Pad



SM324

Carter Braines Support Pad



SM325

Tudor Edwards Support Pad

OPERATING THEATRE ACCESSORIES



SM326

Crawford Wedge

560mm x 280mm x 100mm deep.



SM327

Narrow Arm Table Pad

570mm x 190mm.



SM328

Square Arm Table Pad


405mm x 440mm.



SM329

X-Ray Arm Table Pad

685mm x 480mm



SM330

Abduction Support

For use with a 75mm mattress.

SM331

Abduction Support

For use with a 50mm mattress.

OPERATING THEATRE ACCESSORIES

SM332 **Patient Restraint Strap**
Complete with large retention buckle,
OR Velcro Fixing

SM333 **Patient Restraint Strap**
Complete with stainless steel clamp



SM334 **Anaesthetic Screen Sleeve**
810mm long.



SM335 **Lithotomy Pole Sleeves pair**
460mm long.

SM336 **Lithotomy Pole Sleeve pair**
560mm long.

SM337 **Padded Slings pair**
For Lithotomy poles
with or without billet clips.

SM338 **Canvas Slings pair**
For Lithotomy poles
with or without billet clips.

SM339 **Anti-Static Slings pair**
For Lithotomy poles
with or without billet clips.

OPERATING THEATRE ACCESSORIES



SM340

Montreal Laminectomy Support



SM341

Sacral Support

75mm deep.



SM342

Charnley Type Wedge

SM343

**Leg Trough Covered
Foam Leg Trough Uncovered**

SM344

Ankle/ Arm Protector

SM345

J Shaped Shoulder Rest Pads pair

SM346

Sand Bags

Sizes made to customer's requirements.
To calculate the cost of a standard sand bag
e.g.: £0.21 per square inch 5" x 7".

SAFETY MED LTD

UNIT 1A, HALLAM BUSINESS CENTRE, HALLAM STREET, HEAVILEY, STOCKPORT,
CHESHIRE. SK2 6PT. TEL: 01614773163 FAX: 01614774380

E-mail: sales@safetymedltd.co.uk

www.safetymedical.co.uk